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- Generally, the englobing of the sample is necessary for these reasons:

- The sample is small or it has an irregular shape therefore the handling is difficult.

It is necessary to examine the edge of the sample for the analysis of the thickness of electrochemical coatings, decarburization, heat treated surface, performing microhardness, etc. ...

Facilitate, under certain conditions, the identification of the sample and archiving.

For englobing the sample, depending on the requirements, it is possible to use the cold resin or the hot resin; for the latter, you have to use a mounting press, which, reaching high temperature, and guaranteeing the necessary pressure, hardens the resin making the metallographic sample.

Hot mounting allows to obtain a planar parallel sample. Hot resins are suitable for samples capable of withstanding high pressure and temperatures (~175°)

The cold resin consists in two or more components mixing together (resin + liquid) which are poured into a mold of appropriate size in which the sample is placed. Hardening takes place by chemical reaction with evolution of heat.

Cold resins are suitable for fragile and sensitive samples, which cannot be subjected to high pressure and temperatures.

On the other hand, vacuum cold mounting eliminates any possible air bubbles.

In porous materials, vacuum allows the resin to penetrate into the micro-cavities, compacting the specimen and improving it during the polishing phase.

The Hitech Europe metallographic hot mounting presses line meets the needs of both research laboratories (Universities, Technical Schools, Research Institutes, etc.) and big or small/medium businesses.

The Hitech Europe metallographic hot mounting presses are all fully automatic cycle, operated by hydraulic or pneumatic systems.

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Power 1200 W.
Power supply 220V 50/60Hz 1200 W
Dimensions: 380x450x600 mm
Weight: 42 kg

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Power 2000 W
Power supply 220V 50/60Hz 2000 W
Dimensions: 380x450x600 mm
Weight: 45 kg.

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You can minimize shrinkage and improve edge retention by cooling the mount to room temperature under pressure before removing it from the mounting press

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Undured mounts can be caused by too much moisture in the mounting compound.

Make sure to properly close the container between uses.

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Radial splitting of mounts is often caused by sharp edges on the samples or by samples that are too large for the mold.

Round off sharp corners and move the specimen farther from the edge of the mount.

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Bulging or soft mounts are caused by insufficient cure times. Increase the cure time.

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Unfused or frosted mounting compound is often a sign of insufficient molding temperature or pressure.

Ensure the temperature and pressure settings on the mounting press match the recommendations for the compression compound being using.

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
The presence of cracks on the surface of the sample is due to excessive cooling rate.

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
The presence of 'flakes' inside the specimen is due to too low temperatures or the too short duration of the cycle.

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



SCHEDA TECNICA / TECHNICAL FEATURES

MODEL CODE	EP16 02-22-0A16-C	EP50P 02-22-0P50-C	EP50PT 02-22-PA50-C	EP50T 02-22-0A50-C
Steel frame stainless steel screws	Yes			
Cover	Resin reinforced shockproof / corrosion			Steel
Watt	1200 W	2000 W		
Mould dimensions	1" (25,4mm) 30mm 1 1/4" (31,8mm) 1 1/2" (38,2 mm) 40mm	1" (25,4mm) - 30mm 1 1/4" (31,8mm) 1 1/2" (38,2 mm) - 40mm - 50mm		
Electronics	Thermoregulator		Color Touch Screen	Color Touch Screen
Cycle time	Yes			
Interchangeable Mould	No	Yes		
Closure Mould	Screw	Bayonette		
Power Supply	230 V - AC			
Cooling system	Connect to water source			
Piston driven by	Pneumatic system			Hydraulic motor
Dimensions (mm.)	380x450x600			400x500x650
Weight	42 kg	50 kg	81 kg	
Manual tool set CE Certificate	Yes			

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Tecmet 2000 supplies a wide range of high quality consumables for the preparation of metallographic samples: Abrasive cut off wheels / Diamond blades / Mounting accessories / Hot and cold mounting compounds / Consumables for finishing / Abrasive papers / Self-adhesive diamond discs / Abrasive suspensions / Polishing products / Polishing cloths / Magnetic system

HE HITECH
MILANO EUROPE
MADE IN ITALY BY:

TECNET
QUALITY CONTROL EQUIPMENT 2000

tecmet 2000 srl via Monferrato, 5 - 20094 Corsico(MI) -Italy
tel+39 02 486 00 385 -tel. +39 02 3663 8452
E mail: tecmet@tecmet2000.it www.tecmet2000.it